

# Detection of 11-Nor-9-carboxy- $\Delta^9$ -Tetrahydrocannabinol in Urine using Thomson eXtreme|FV<sup>®</sup>

Taylor Yurasits  
*York College*

Crystal Xander BS  
*Health Network Labs*

Dear Fritch PhD  
*Health Network Labs*

Follow this and additional works at: <http://scholarlyworks.lvhn.org/research-scholars-posters>

---

## Published In/Presented At

Yurasits, T., Xander, C., Fritch, D., (2017, July, 31) *Detection of 11-Nor-9-carboxy- $\Delta^9$ -Tetrahydrocannabinol in Urine using Thomson eXtreme|FV<sup>®</sup>*. Poster presented at LVHN Research Scholar Program Poster Session, Lehigh Valley Health Network, Allentown, PA.

This Poster is brought to you for free and open access by LVHN Scholarly Works. It has been accepted for inclusion in LVHN Scholarly Works by an authorized administrator. For more information, please contact [LibraryServices@lvhn.org](mailto:LibraryServices@lvhn.org).

# Detection of 11-Nor-9-carboxy- $\Delta^9$ -Tetrahydrocannabinol in Urine using Thomson eXtreme|FV®

Taylor Yurasits\*<sup>1</sup>, Crystal Xander B.S.<sup>2</sup> and Dean Fritch Ph.D.<sup>2</sup> York College of Pennsylvania, York, PA <sup>2</sup>Health Network Laboratories, Allentown, PA

Lehigh Valley Health Network, Allentown, Pennsylvania

## BACKGROUND / INTRODUCTION

- To determine and validate a more time effective extraction method for quantifying 11-Nor-9-carboxy- $\Delta^9$ -Tetrahydrocannabinol (THC-COOH) in urine
- Tetrahydrocannabinol or THC is one of the most commonly used illicit drugs, but has some medical acceptance
- The main metabolite of THC in urine is THC-COOH in a glucuronide form
- Solid Phase Extraction (SPE) using the Zymark Rapid Trace extracts THC-COOH in urine from other compounds by its chemical and physical properties.
- The Thomson eXtreme|Filter Vial® (TV) filters out impurities from the sample and can run instantly on the LC/MS/MS.
- Thomson eXtreme|FV® 0.2  $\mu$ m PVDF was used to lower sample preparation time.

## METHODS

### Quality Control

- The controls and patient samples were run on a 5 point calibration curve for validation. The curves were spiked at 15, 30, 75, 750 and 3000 ng/mL. The standards were obtained from Cerilliant, Round Rock, Texas, at concentrations of 100  $\mu$ g/mL or 1.0 mg/mL.

### Sample Preparation

- Urine samples were hydrolyzed with 90  $\mu$ l of  $\beta$ -glucuronidase and incubated at 55°C for 30 minutes. 200  $\mu$ l of sample and 200  $\mu$ l of 2% methanol in HPLC grade water were placed into Thomson eXtreme|FV® 0.2  $\mu$ m PVDF. Samples were pressed to remove impurities before LC/MS/MS analysis.


### Instrumentation

- All hydrolyzed samples were run on a Shimadzu Prominence Sciex 5500 Qtrap LC/MS/MS.



Figure 1: Thomson eXtreme|FV® during filtering process.

## OUTCOMES




Graph 1: A calibration curve for the analyte, THC-COOH, over 5 days had an R<sup>2</sup> value of 0.99592.

| Between Run    |               |                |
|----------------|---------------|----------------|
| Sample         | TV % Accuracy | SPE % Accuracy |
| LOQ            | 103.3         | 103.3          |
| LOW CONTROL    | 92.2          | 95.4           |
| MIDDLE CONTROL | 91.1          | -              |
| HIGH CONTROL   | 84.6          | 84             |
| LOW STABILITY  | 89.6          | -              |
| HIGH STABILITY | 88.5          | -              |

Figure 2: Accuracy of quality control samples between all 5 days.

### Time Evaluation Per Sample



Graph 2: Comparison of time between two extraction methods.

## RESULTS

- Figure 1 depicts the process of the Thomson eXtreme|FV® removing macroparticles from the sample that are not necessary for sample injection.
- Based off of Graph 1, the method reached accurate and precise results with an R<sup>2</sup> value of 0.99592 over 5 days.
- Figure 2 shows the limit of quantification (LOQ) for both methods have the same accuracy.
- The low and high control are very close in precision between both methods.
- Graph 2 shows that the new Thomson vial method takes 5 minutes per samples and the old Rapid Trace method took 9 minutes per sample.

## CONCLUSIONS

- The accuracy and precision did not change the LOQ when using the new Thomson eXtreme|FV® 0.2  $\mu$ m PVDF method compared to the old SPE Rapid Trace method. The high and low controls accuracy was only a minor difference, but still having correctness.
- By using the Thomson eXtreme|FV® method, the time decreases by 4 minutes. Increasing the extraction efficiency increases the turn around time for reporting of patient samples.

### Acknowledgments

I would like to thank Crystal Xander B.S. for her support and guidance throughout the project and Dean Fritch Ph.D. for assisting in the data analysis.

### References

- <sup>1</sup>"Cerilliant-Certified Reference Materials | Certified Reference Standards | Analytical Standards". *Cerilliant.com*. 2017.
- <sup>2</sup>*IMCSzyme*. 2017.
- <sup>3</sup>"Filter Vials | Thomson Instrument Company". *Htslabs.com*.
- <sup>4</sup>Science, Live. "What Is THC?". *Live Science*. 2017

© 2017 Lehigh Valley Health Network

610-402-CARE LVHN.org

