

Setting Optimal Positive Expiratory Pressure Utilizing Transpulmonary Monitoring and a Pressure Volume Measurement.

Kenneth Miller MEd, RRT-NPS

Lehigh Valley Health Network, Kenneth.Miller@lvhn.org

Matthew Reis RRT, ACCS

Lehigh Valley Health Network, Matthew_J.Reis@lvhn.org

Dave Marth RRT

Lehigh Valley Health Network, David.Marth@lvhn.org

Nirupama Kakumanu MD, FCCP

Lehigh Valley Health Network, Nirupama.Kakumanu@lvhn.org

Follow this and additional works at: <https://scholarlyworks.lvhn.org/patient-care-services-nursing>



Part of the [Nursing Commons](#), and the [Pulmonology Commons](#)

Published In/Presented At

Miller, K. Reis, M. Marth, D. Kakumanu, N. (2017, October 3). *Setting Optimal Positive Expiratory Pressure Utilizing Transpulmonary Monitoring and a Pressure Volume Measurement*. Poster Presented at: The American Association of Respiratory Care (AARC), Innapolis, Indiana.

This Poster is brought to you for free and open access by LVHN Scholarly Works. It has been accepted for inclusion in LVHN Scholarly Works by an authorized administrator. For more information, please contact LibraryServices@lvhn.org.

Setting Optimal Positive Expiratory Pressure Utilizing Transpulmonary Monitoring and a Pressure Volume Measurement


Kenneth Miller MSRT, MEd, RRT-ACCS, FAARC, Matthew Reis RRT-ACCS, David Marth RRT, Nirupama Kakumanu MD, FCCP
Respiratory Care Department, Lehigh Valley Health Network, Allentown, Pennsylvania

INTRODUCTION

- Setting optimal positive end expiratory pressure (PEEP) is beneficial in improving gas exchange, promoting alveolar inflation, and minimizing ventilator induced trauma.
- The setting of PEEP has historically been done many different ways, FIO₂/PEEP algorithmic approach, compliance assessment, etc.
- Recently, we have employed the use of trans-pulmonary monitoring along with a pressure/volume measurement (P/V) to set the optimal PEEP.
- By utilizing both of these strategies the optimal PEEP can be set to maximize clinical objectives.

METHODS

- We inserted an esophageal balloon (Cooper Surgical, USA) to monitoring TranspE.
- We performed a pressure/volume(G-5 Hamilton Medical, Switzerland) to determine the lower inflection point on twelve patients diagnosed with severe to moderate ARDS. (Berlin definition)
- All patients were sedated and placed on a paralytic infusion.
- Esophageal balloon placement was verified via ventilator waveform and chest X-ray.
- Our goal was to set the PEEP to maintain a TranspE between -2 to +2cm/H2O and 2 cmH2O above the lower inflection point via the P/V measurement (G-5 P/V Tool).



Post esophageal balloon placement.

Peak airway pressure minus esophageal pressure equals:
Transpulmonary pressure (pressure the lung is receiving)
Goal <25 cm

Note peak airway pressure is 42cm
Esophageal pressure 21cm thus transpulmonary pressure
(the lung receiving) only 19cm

PEEP set at 14cm
Based on PtransE
Between =2 to +2cm



RESULTS

- In nine of the patients the optimal PEEP was the same via the TranspE and P/V measurements. (Table 1)
- In three patients, all >125kg there was a greater than 4cm/H2O difference between the TranspE and P/V measurement lower inflection (LIP) point regarding the optimal PEEP setting.
- On these three patients we set the PEEP between the two measurements and assessed clinical performance.
- Post PEEP setting in all twelve patients Transpl was maintained <25cm/H2O and hemodynamic performance remained stable.

Table 1. TranseP and P/V Measurements

Wt/Kg	TranspE cm/H ₂ O	LIP cm/H ₂ O	PEEP Set cm/H ₂ O
95	-1.5	14	16
150	+6	20	17
80	+1	12	14
55	-2	12	14
75	+2	10	12
110	-2	16	18
175	-6	16	22
130	-4	18	20
92	+1	12	14
89	-2	12	14
105	+1.5	14	16
72	-1	10	12

CONCLUSIONS

- Based on our experience, we believe utilizing both transpulmonary monitoring and pressure volume measurements can determine the optimal PEEP.
- Factors in chest wall impedance may cause a difference in these measurements and enhanced clinical assessment should be utilized to determine the optimal PEEP in that patient population.
- The effect of these measurements on mortality and morbidity still needs to be determined.

© 2017 Lehigh Valley Health Network

610-402-CARE LVHN.org